


UNICORD[®] Endoscopy Cable Ordering Guide

**PRODUCTS
FOR
MEDICINE**

ORDER NUMBER




Example Order Number


AC-2010-ST

2.0 mm, 10 ft cable with ACMI Light Source Connector and Storz Instrument Connector

2. Specify Bundle Diameter

20	2.0 mm	 <p>To Determine: Measure fiber diameter at end of cable.</p>
35	3.5 mm	
50	5.0 mm	

3. Specify Cable Length

60	6 Feet	 <p>Cable Length is measured from fiber bundle terminations.</p>
75	7.5 Feet	
10	10 Feet	

4. Specify Instrument Connector

(800 Series)

AC	ACMI, STRYKER, ZIMMER, EDER, LUXTEC, ASPEN, AFI, ZEISS USA, WEHMER, KLI, V. MUELLER, CODMAN, CONCEPT
PI	PILLING
WO	WOLF, DYONICS, WECK, LINVATEC
ST	STORZ
DV	DESIGNS FOR VISION
OL	OLYMPUS

1. Specify Light Source Connector

(900 Series)

AC	ACMI, STRYKER, ZIMMER, EDER, LUXTEC, ASPEN, AFI, ZEISS USA, WEHMER, KLI, V. MUELLER, CODMAN, CONCEPT
PI	PILLING
WO	WOLF, DYONICS, WECK, LINVATEC
ST	STORZ
DV	DESIGNS FOR VISION
OL	OLYMPUS

Note:

All instrument connectors utilize snap-ring design. Alternate cable lengths available upon request. Please call for further information.

Product and/or company names are trademarks or registered trademarks of those respective companies.

**PRODUCTS
FOR
MEDICINE**

1201 East Ball Road Suite H, Anaheim, CA 92805

<http://www.productsformedicine.com>

Toll Free: (800) 333-3087

Fax: (714) 991-8200